



# Welcome to our hands-on tour of the single cell analysis journey

From Cloud to Loupe to community-developed tools

---

CG000586 RevB

# Goals and objectives

## Primary goal

- **Demystify the data analysis process with a hands-on, active learning approach.**

## Target audience


- Those who have or will soon have 10X data and have limited to no bioinformatics experience

## Secondary goals

### **TL;DR: Help you move forward in analysis**

- Help you see a path for you through the analysis journey
- Bring awareness to useful analysis resources, such as software documentation and analysis guides
- Get you started in raw data processing using 10x cloud analysis
- Show the expanded functionality of our Loupe Browser
- Introduce the amazing ecosystem of community developed tools for analyzing 10X data

# Analysis Journey: Orientation



# Agenda

## Raw data processing and QA

- Introduction
- 10x Genomics Cloud Analysis
- Overview of scRNA-seq Data Analysis
- Quality Assessment

## Downstream analysis and visualization

- Loupe Analysis
- RNA Velocity: Community Developed Tools
- Planning your Individual Journey



# Introduction to the neutrophil dataset: analysis begins at design

---

Hands on tour of the single cell analysis journey

# Outline

- Analysis begins at design
- The data we are using in this series
- Purpose of the experiment
- Neutrophil challenges
- Analysis plan

## Analysis begins at design

- The design enables analysis
- What samples, conditions, and data do you need?
- The planning call is the cheapest part of the experiment



<https://imgflip.com/memegenerator/One-Does-Not-Simply>

# The data we are using today

10x Genomics Support / Single Cell Gene Expression / Documentation / Sample Prep /

## Neutrophil Analysis in 10x Genomics Single Cell Gene Expression Assays



Technical Note, CG000444

CG000444\_Neutrophil Analysis\_TN\_RevA .pdf

View and download

### Document Type

Technical Note

### Last Modified

October 25, 2021

Neutrophils are the most abundant white blood cells in circulation. They are the first responders of the innate immune system that release cytotoxic compounds to kill bacteria and phagocytose foreign particles. Analysis of neutrophils and granulocytes is challenging in single cell transcriptomic data due to the low RNA content in these cells along with relatively high levels of RNases and other inhibitory compounds. Furthermore, neutrophils are sensitive to degradation after collection, requiring careful sample handling for preservation. This Technical Note describes the successful isolation and analysis of neutrophil populations from whole blood.

<https://www.10xgenomics.com/support/single-cell-gene-expression/documentation/steps/sample-prep/neutrophil-analysis-in-10-x-genomics-single-cell-gene-expression-assays>




# Purpose of the experiment

- Show neutrophils are detected by the assay
- Show that the biology is preserved




# Neutrophil challenges

- Neutrophils are hard to collect on their own
- Neutrophils have lower RNA content and express fewer genes
- Neutrophils have elevated levels of intron retention



# Neutrophil analysis plan

Our journey through analysis



## Next up:

10x Genomics Cloud Analysis: from FASTQs to Quantitative Biological Information with Cell Ranger

### Topics

- Analysis begins with design: the goals of the experiment will determine how the data should be processed and analyzed.
- Introduction to the neutrophil data set

### Next

- Signing in and creating a project
- Using the web based FASTQ uploader
- Creating a new analysis
- Downloading output files