

# Welcome to our hands-on tour of the single cell analysis journey

from Cloud to Loupe to community-developed tools

---

7/20/2022

University of Michigan, Ann Arbor

CG000586 RevA

# Welcome! And thank you to our Hosts!

- We are excited to be here and to meet you
- Special thanks to Olivia (Liv) Koues and Chris Gates



**Genomics Core**  
734-764-8531  
okoues@umich.edu

Rm 119, 2800 Plymouth Road  
NCRC B14  
Ann Arbor MI 48103-2800



**Bioinformatics Core**  
734-998-9249  
cgates@umich.edu

North Campus Research Complex  
NCRC B22  
Ann Arbor, MI 48109

# Agenda

Morning: Raw data processing and QA

- Introduction and housekeeping
- 10x Genomics Cloud Analysis
- Introduction to scRNA-seq Data Analysis
- Break
- Quality assessment

Lunch

Afternoon: Downstream analysis and visualization

- Loupe analysis
- Break
- RNA velocity: third party tools
- Planning your individual journey
- Wrap-up

# Code of conduct


- Ask questions
- Be helpful
- Call out if you get stuck
- Reminder: pre-workshop survey



# Workshop team and helpers



Dann



Juan



Lisa



Matt



Mike



Nur



Ross

# Next up:

10x Genomics Cloud Analysis: from FASTQs to Quantitative Biological Information with Cell Ranger

<https://web-frontend-6681okbdp-10x-genomics.vercel.app/resources/analysis-guides/workshop-cloud-walkthrough>

## Topics

- Analysis begins with design: the goals of the experiment will determine how the data should be processed and analyzed.
- Signing in and creating a project
- Using the web based FASTQ uploader
- Creating a new analysis
- Downloading output files



# Analysis Journey: Orientation

## Experimental Design

Analysis begins with design. Choose assay types and define controls & replicates to set your course.

## Lab Workflow

Use analysis-aware execution to avoid road-blocking technical artifacts.

## Raw Data Processing


Transform your sequencing data into quantitative biological information.

## Quality Control

Assess the quality of your data at multiple stages of analysis.

## Visualization & Interpretation

Arrive at biological insights and address research questions with Loupe Browser and community-developed tools.



# Introduction to the neutrophil dataset: analysis begins at design

---

7/20/2022

Hands on tour of the single cell analysis journey



# Outline

- Analysis begins at design
- The data we are using today
- Purpose of the experiment
- Neutrophil challenges
- Analysis plan

# Analysis begins at design

- The design enables analysis
- What samples, conditions, and data do you need?
- The planning call is the cheapest part of the experiment



# The data we are using today

10x Genomics Support / Single Cell Gene Expression / Documentation / Sample Prep /

## Neutrophil Analysis in 10x Genomics Single Cell Gene Expression Assays



Technical Note, CG000444

CG000444\_Neutrophil Analysis\_TN\_RevA .pdf

View and download

### Document Type

Technical Note

### Last Modified

October 25, 2021

Neutrophils are the most abundant white blood cells in circulation. They are the first responders of the innate immune system that release cytotoxic compounds to kill bacteria and phagocytose foreign particles. Analysis of neutrophils and granulocytes is challenging in single cell transcriptomic data due to the low RNA content in these cells along with relatively high levels of RNases and other inhibitory compounds. Furthermore, neutrophils are sensitive to degradation after collection, requiring careful sample handling for preservation. This Technical Note describes the successful isolation and analysis of neutrophil populations from whole blood.

<https://www.10xgenomics.com/support/single-cell-gene-expression/documentation/steps/sample-prep/neutrophil-analysis-in-10-x-genomics-single-cell-gene-expression-assays>


# Purpose of the experiment

- Show neutrophils are detected by the assay
- Show that the biology is preserved



# Neutrophil challenges

- Neutrophils are hard to collect on their own
- Neutrophils have lower RNA content and express fewer genes
- Neutrophils have elevated levels of intron retention



# Neutrophil analysis plan

Our journey through analysis

